



MERCURIUS™

**High sensitivity low-input  
BRB-seq Service**

---

**Sample Submission Guidelines**

January 2025

## Sample submission guidelines at a glance

1. Transfer the RNA samples to a 96- or 384-well PCR plate following the instructions below and store it at -80°C before shipment.
2. Fill in the Sample Submission Form (**SSF**) and **check all the boxes** in the Sample submission checklist below; send both files to **orders@alitheagenomics.com**.

Please be aware that any inconsistency may lead to delays or additional fees.

3. Ship the samples on dry ice and send us the shipment tracking number.

## Sample submission checklist

- Sample Submission Form (SSF)** is filled properly with a **unique sample ID**. Consider adding a suffix for technical replicates (e.g., XX\_rep1, XX\_rep2, etc.). Make sure that information about **all the samples shipped is provided in the SSF**.
- The minimum sample volume is **15 µL** per well.
- The concentration of all samples is uniform, and the A260/230 ratio is assessed and provided in the SSF.
- The **minimum number** of samples in each group (to be pooled together) is **16**.
- Samples are provided in the 96- or 384-well RNase/DNase-free PCR plates. **Samples in tubes cannot be processed**.
- Plates with samples are clearly labeled with the same Plate ID as in the sample submission form.
- Plates are well sealed with adhesive and temperature-resistant seal (aluminum is ideal).

## Essential considerations for input material

### Samples quantity and integrity

- The tested range of total RNA amount is 100 pg – 500 ng (per well).
- The minimum sample volume should be 15 µL.
- The sample concentration should be at least 2 ng/µL (for the QC procedure).
- The recommended RIN number is > 6.

### Samples purity

- RNA samples extracted with TRIzol, phenol, chloroform, or guanidine are prone to residual contamination with organic solvents, which considerably decreases cDNA yield.
- To ensure the high purity of RNA, assess the 260/230 ratio for all (or at least for 10% of samples) from the same RNA purification batch.
- The 260/230 ratio values should be between 1.8 and 2.2.

### Samples uniformity

- To ensure an even distribution of reads after sequencing, the RNA amount, integrity, and 260/230 values of the starting RNA samples must be as uniform as possible, with a max 10% variation.
- To obtain such uniform amounts, we recommend the following:
  - Use dye-based methods for RNA quantification (e.g., Qubit Quant-iT or RiboGreen).
  - Dilute samples to obtain the same RNA concentration in all wells ( $\pm 10\%$ ).
  - Re-measure the RNA concentration of all samples to confirm uniform concentration.
  - Ensure that the 260/230 ratio is between 1.8 and 2.2 and the RIN value is > 6 across the samples.

### Batch-effect and sample replicates

- The RNA extraction protocol can produce considerable technical variation across the samples; therefore, it is strongly recommended that RNA extraction be performed in a single batch.
- If the differential expression (DE) analysis is planned, the respective RNA samples should be included in the same library. Comparing samples from different libraries can be biased.
- It is highly recommended that at least 3 (or more) biological replicates be included for the reliability of the experimental setup.

### Samples preparation

1. Label a new 96- or 384-well RNase/DNase-free PCR plate.
2. Pipette the RNA samples to the new 96- or 384-well PCR plate according to the filled Sample Submission Form. Follow the column-based direction (column 1, then column 2, etc.).
3. Seal the 96- or 384-well PCR plate with an aluminum seal and briefly spin it down
4. Store the samples at  $-80^{\circ}\text{C}$  before shipment.

## Alithea Genomics SA

**Phone** +41 78 830 31 39

**Email** [info@alitheagenomics.com](mailto:info@alitheagenomics.com)

**Web** [www.alitheagenomics.com](http://www.alitheagenomics.com)

**Address**

Route de la Corniche, 8  
1066 Épalinges  
VD, Switzerland