

ALITHEA
GENOMICS

MERCURIUS™

**Organoid DRUG-seq
Library Preparation kit
for 96, 384, and 1'536 samples**

PN 10861, 10871, 11061, 11071

User Guide

v.0.9.1.revA December 2024

FOR RESEARCH USE ONLY | ALITHEAGENOMICS.COM

Related Products

Kit name	Kit PN	Modules	Module PN
Mercurius Organoid DRUG-seq Library Preparation 96 kit	10861	Barcoded Oligo-dT Adapters Module 96 samples	10513
		Organoid DRUG-seq Library Preparation and UDI Module 96 samples	10545
Mercurius Organoid DRUG-seq Library Preparation 384 kit	10871	Barcoded Oligo-dT Adapters Module 384 samples	10515
		Organoid DRUG-seq Library Preparation and UDI Module 384 samples	10571
Mercurius Organoid DRUG-seq Library Preparation 384 (4x 96) kit	11061	Barcoded Oligo-dT Adapters Module 4x 96 samples	10513
		Organoid DRUG-seq Library Preparation and UDI Module 4x 96 samples	10565
Mercurius Organoid DRUG-seq Library Preparation 1'536 (4x 384) kit	11071	Barcoded Oligo-dT Adapters Module 4x 384 samples	10515
		Organoid DRUG-seq Library Preparation and UDI Module 4x 384 samples	10566

Table of Contents

TABLE OF CONTENTS	2
KIT COMPONENTS	3
Reagents supplied	3
Additional required reagents and equipment (supplied by the user).....	4
PROTOCOL OVERVIEW AND TIMING	5
PROTOCOL WORKFLOW	6
PART 1. PREPARATION OF ORGANOID LYSATE SAMPLES	7
1.1. Essential considerations for input material	7
1.2. ERCC Spike-in Controls (Optional).....	7
1.3. Organoid pellet preparation	7
1.4. Organoid lysate preparation.....	8
PART 2. LIBRARY PREPARATION PROTOCOL	9
2.1. Reverse transcription	9
2.2. Sample pooling and column purification	10
2.3. Free primer digestion	12
2.4. Second strand synthesis.....	13
2.5. Tagmentation	14
2.6. Library indexing and amplification.....	15
2.7. Library quality control	17
PART 3. LIBRARY SEQUENCING	19
3.1. Sequencing with Illumina instruments.....	19
3.2. Sequencing with Element AVITI instrument.....	20
PART 4. SEQUENCING DATA PROCESSING	21
4.1. Required software	21
4.2. Data processing	21
4.2.1. Merging <i>fastq</i> files from individual lanes and/or libraries (Optional).....	21
4.2.2. Sequencing data quality check.....	22
4.2.3. Preparing the reference genome.....	22
4.2.4. Aligning to the reference genome and generation of count matrices	23
4.2.5. Generating the count matrix from .mtx file	24
4.2.6. Generating the read count matrix with per-sample stats (Optional)	24
4.2.7. Demultiplexing bam files (Optional).....	24
APPENDIX 1. ERCC SPIKE-IN CONTROL	26
Organoid Lysis Buffer (OLB) preparation with ERCC	26
APPENDIX 2. COMPATIBLE ILLUMINA INSTRUMENTS	27
CONTACT INFORMATION	28

Kit Components

Reagents supplied

Barcoded Oligo-dT Adapters Set V5 Module

Component Name	Label	Amount provided per kit				Storage
		96 samples (PN 10861)	384 samples (PN 10871)	4x 96 samples (PN 11061)	4x 384 samples (PN 11071)	
Plate with 96 barcoded oligo-dT primers, set V5D (PN 10419)	96 V5D OdT	1 plate	-	4 plates	-	-20°C
Plate with 384 barcoded oligo-dT primers, set V5D (PN 10420)	384 V5D OdT	-	1 plate	-	4 plates	-20°C
Aluminium Seal	-	3 pcs	3 pcs	12 pcs	12 pcs	-20°C/RT

Organoid DRUG-seq Library Preparation and UDI Module

Component Name	Label	Cap colour	Volume (µL)				Storage
			96 samples (PN 10545)	384 samples (PN 10571)	4x 96 samples (PN 10565)	4x 384 samples (PN 10566)	
RT Enzyme	RTE	magenta	11	43	43	4x 43	-20°C
RT Buffer	RTB	magenta	550	1100	2x 1100	4x 1100	-20°C
Exonuclease I Enzyme	EXO	purple	10	10	10	10	-20°C
Exonuclease Buffer	EXB	purple	20	20	20	20	-20°C
Second Strand Enzyme	SSE	orange	20	20	20	20	-20°C
Second Strand Buffer	SSB	orange	45	45	45	45	-20°C
Tagmentation Enzyme	TE	red	12	12	12	12	-20°C
Tagmentation Buffer	TAB	red	40	40	40	40	-20°C
Library Amplification Mix	LAM	green	200	200	200	200	-20°C
UDI Adapter Mix 1	MQ.UDI.1	transparent	10	10	10	10	-20°C
UDI Adapter Mix 2	MQ.UDI.2	transparent	10	10	10	10	-20°C
UDI Adapter Mix 3	MQ.UDI.3	transparent	10	10	10	10	-20°C
UDI Adapter Mix 4	MQ.UDI.4	transparent	10	10	10	10	-20°C
Dissociation Reagent	DIR	yellow	1575	2x 1575	4x 1575	8x 1575	-20°C
Organoid Lysis Buffer	OLB	yellow	1050	2x 1050	4x 1050	8x 1050	-20°C
RNase Inhibitor	OIN	yellow	263	2x 263	4x 263	8x 263	-20°C

Additional required reagents and equipment (supplied by the user)

Plasticware	Manufacturer	Product number
15 mL conical tubes	Greiner	188271
0.2 mL 8-Strip Non-Flex PCR Tubes	Starlab	I1402-3700
Zymo-Spin IICG Columns (optional)	Zymo	C1006-50-G
25 mL reservoir for spin columns	Zymo	C1039-25
Disposable Pipetting Reservoir 25mL polystyrene	Integra	4382
Solution reservoir 150 mL polystyrene (optional)	Integra	6318

Reagents	Manufacturer	Product number
DNA Clean and Concentrator-5 kit	Zymo	D4014
SPRI AMPure Beads	Beckman Coulter	A63881
or		
CleanNA Beads	CleanNA	CNGS0050D
Qubit™ dsDNA HS Assay Kit	Invitrogen	Q32851
High Sensitivity NGS Fragment Analysis Kit	Agilent	DNF-474
SYBR Green	Thermofisher	S7563
Ethanol, 200 proof	-	-
Nuclease-free water	Thermo Fisher	A57775
DPBS, Cell culture grade	Gibco	10010023
ERCC RNA Spike-In Mix	Thermo Fisher	4456740

Equipment	Manufacturer	Product number
Benchtop centrifuge for plates	-	-
Benchtop centrifuge for 1.5 mL tubes	-	-
Single and Multichannel pipettes	-	-
Fragment Analyser / Bioanalyzer / TapeStation	Agilent	M5310AA
Qubit™	Invitrogen	Q33238
Magnetic stand for 0.2 mL tubes	Permagen	MSR812
Magnetic stand for 1.5 mL tubes	Permagen	MSR06
12-channel pipette, 0.5-12.5 µL VIAFLO or similar	Integra	4631
12-channel pipette, 5-125 µL VIAFLO or similar	Integra	4632
8-channel adjustable tip spacing pipette, VOYAGER, 2 – 50 µL	Integra	4726
Pipetboy	Integra	155 000
Vacuum manifold Qiagen Qiavac 24 Plus or similar	Qiagen	19413
Vacuum pump	-	-
Real-Time PCR Instrument (optional)	-	-
VIAFLO instrument (optional)	Integra	6001
VIAFLO 96/384 channel pipetting head, 0.5-12.5 µL (optional)	Integra	6101/6131

Protocol Overview and Timing

The MERCURIUS™ Organoid DRUG-seq kits allow the preparation of Illumina-compatible 3' RNA sequencing libraries for up to 1'536 organoid samples in a time and cost-efficient manner. The kits include a three-component lysis buffer for efficient organoid and spheroid digestion and cell lysis. Obtained crude lysates can be used directly in the RT reaction, skipping the tedious RNA purification step.

The kits are provided in the following formats:

Kit format	PN	PCR plate format	Maximum number of samples in one pool	Maximum number of samples processable	Number of UDI Libraries
96-sample	10861	96WP	96	96	4
384-sample	10871	384WP	384	384	4
4x 96-sample	11061	96WP	96	384	4
4x 384-sample	11071	384WP	384	1'536	4

Every kit contains barcoded MERCURIUS™ Oligo-dT primers designed to tag polyA+ RNA from lysate samples with individual barcodes during the first strand synthesis reaction, allowing the pooling of the resulting cDNA samples from each experimental group into a single tube for streamlined sequencing library preparation.


The BRB-seq technology can be used to generate high-quality sequencing data from 2'000 – 50'000 mammalian cells per well. Notably, the kit can be used to pool any number of samples up to the capacity of the provided plate (96 or 384) with two considerations:

- The total cell number per pool should be at least 80'000 cells.
- Pooling less than eight samples may result in low-complexity reads during sequencing, decreasing the overall sequencing quality. If necessary, the latter can be improved by increasing the proportion of PhiX spike-in control during sequencing (see [Part 3](#)).

Each library indexing is performed using a Unique Dual Indexing (UDI) strategy, which minimizes the risk of barcode misassignment after NGS. The kits contain four UDI adapters. Every adapter can be used to prepare an individual library. Libraries with different UDI adapters can be pooled and sequenced in a single flow cell.

[Figure 1](#) provides an estimation of the time required to accomplish each step of the protocol.

Protocol Workflow



Overall time

- Incubation time: 3h30-3h45.
- Hands-on time: 2h20-3h10.
- QC time: 5min-1h20 (depending on the instrument used: Qubit or Fragment Analyzer).
- Sequencing time: 19-48 hours (depending on instrument and sequencing run settings).

* The first provided time is for a 96 well-plate, followed by the time for a 384 well-plate.

Figure 1 Schematic illustration of the protocol workflow

Part 1. PREPARATION OF ORGANOID LYSATE SAMPLES

NOTE: Before starting every step, briefly spin down the tubes and plates before opening them to ensure that all liquid or particles are collected at the bottom of the tube/plate.

1.1. Essential considerations for input material

- The organoid lysate protocol was validated with organoids and spheroids cultured in flat-bottom 96-well plates with Matrigel®, or microtissue-specific 96-well plates (e.g., Akura™, BIOFLOAT™, etc.).
- The recommended input range of cells in organoids or spheroids is 5'000-50'000 cells/well of **96WP** and 2'000-10'000 cells/well of **384WP**.
- Organoids must be plated a few days in advance for the best results.
- Depending on the type of organoids or spheroids (human, mouse, metastatic, or primary cells) and experimental design (e.g., drug treatment, induction of apoptosis, cell cycle arrest, etc.), consider the organoid-specific cell doubling time and potential effect of the treatment on cell quality and quantity.
- The lysis protocol requires vigorous pipetting, specifically if organoids were grown in solidified Matrigel®.
- To ensure an even distribution of reads after sequencing, the amount of starting material must be as uniform as possible. For this, we suggest using automatized cell seeding instruments or double-verified cell counts.
- The organoids lysate used in the RT reaction should be at most 50% of the final volume.

1.2. ERCC Spike-in Controls (Optional)

To ensure the evaluation of the sequencing reads uniformity across the samples and to assess the impact of sample and library preparation steps on it, we recommend the addition of External RNA Controls Consortium (ERCC) Spike-Ins to the lysate buffer (Thermo Fisher, 4456740). Please follow [Appendix 1](#) for detailed information before the lysis step.

1.3. Organoid pellet preparation

At this step, plated organoids or spheroids are washed with DPBS and frozen at -80°C for at least 5 minutes. If possible, snap-freeze the plate with dry ice or liquid nitrogen beforehand.

Preferably perform the procedure of the lysis step just before doing the reverse transcription reaction.

NOTE: The freezing step is required to achieve a higher exon mapping.

- 1.3.1. Gently aspirate culture media from the plate and wash organoids by adding DPBS (per well):
 - **96WP** (Flat-bottom plates with Matrigel®): 80-100 µL DPBS
 - **96WP** (Microtissue-specific plates) and **384WP**: 20µL DPBS
- 1.3.2. Centrifuge at 300x g for 3 min if needed.
- 1.3.3. Gently tap the plate and aspirate as much DPBS as possible without disturbing the structure of the organoids.
- 1.3.4. Seal the plate well with an Aluseal and immediately transfer it to a -80°C freezer for storage. If possible, snap-freeze the plate with dry ice beforehand.

NOTE: If several plates must be processed, perform the procedure with each plate individually one by one to avoid keeping plates at room temperature for a prolonged time.

1.4. Organoid lysate preparation

At this step, frozen organoids are lysed directly in a 96-well plate by adding 1x Cell Lysis Buffer to the wells. The lysates can be used directly for the mRNA selection.

Preparation

Depending on the type of plates, we recommend using:

- **96WP** (Flat-bottom plates with Matrigel®): 30 μ L
- **96WP** (Microtissue-specific plates) and **384WP**: 10 μ L

Prepare a working solution of 1x Organoid Lysis Buffer according to the type and number of wells:

Reagent	96WP (Flat-bottom), μ L		96WP (Microtissue-specific), 384WP, μ L	
	Per well	96 wells +10%	Per well	384 wells +10%
DIR	15	1575	5	2200
OLB	10	1050	3.4	1496
OIN	2.5	263	0.8	352
Water	2.5	262	0.8	352
TOTAL	30	3150	10	4400

Pipette a prepared mix gently a few times, and briefly spin the tube. Keep the mix on ice until further use.

Procedure for organoid lysis

- 1.4.1. Depending on the type of plates, use a multi dispenser to distribute 10 μ L or 30 μ L of a prepared OLB per well.
- 1.4.2. Centrifuge the plate at 300x g for 1 minute to ensure that OLB is uniformly distributed on the surface of each well.
- 1.4.3. Gently pipette the OLB 10-20 times and incubate a plate at +37°C for 15 min. To reduce bubbles, preferably use a multichannel pipette and avoid intensive pipetting.
- 1.4.4. Optionally, perform a second round of pipetting after 5-10 min of incubation.
- 1.4.5. Centrifuge a plate at 600x g for 3 min.
- 1.4.6. Label a new 96- or 384-well PCR plate.
- 1.4.7. Transfer the entire volume of organoid lysates from every well to a corresponding well of the 96- or 384-well PCR plate. Preferably use a multichannel pipette and avoid transferring any cell pellet.
 - Pipette **directly into the plate containing oligos** (follow step 2.1.1); or
 - Pipette **to a new PCR plate** for a future experiment. Ensure proper labelling for easy identification during subsequent use.
- 1.4.8. Seal the 96-well or 384-well PCR plate with an aluminium seal provided and briefly spin it down.
- 1.4.9. The lysates can be used directly for library preparation (see step 2.1, preferably) or safely stored at -80°C.

Part 2. LIBRARY PREPARATION PROTOCOL

2.1. Reverse transcription

Each individual organoid lysate sample is reverse-transcribed at this step using the barcoded oligo-dT primers provided in a 96-well (**96WP**) or a 384-well (**384WP**) plate format depending on the kit type. Subsequently, all barcoded samples can be pooled in one tube.

NOTE: Barcoded oligo-dT primers are provided lyophilized with the addition of dye. The dye has no impact on the enzymatic reactions and is used solely for better visualization of reaction preparation and pooling.

Despite variations in appearance caused by the drying process, wells may exhibit traces of dried dye ranging from dispersed to solid dots on the bottom. The following addition of RT reagents will enable the visualization of red color, confirming the presence of the oligos in all wells.

Preparation

- Thaw the organoid lysate samples on ice.
- Thaw the **RTB** reagent at room temperature and mix well before use.
- Briefly spin down the **96WP** or **384WP** plate containing dried oligo-dT primers. This plate will be referred to as the RT plate.
- Prepare program **1_RT** on the thermocycler (set the lid at 90°C):

Step	Temperature, °C	Time
Incubation	50	30 min
Inactivation	85	10 min
Keep	4	pause

NOTE: All the manipulations with organoid lysates and RT enzyme should be performed in an RNase-free environment, with the use of RNase-free consumables and filter tips, on ice, and using gloves.

Procedure for **96WP** and **384WP**

NOTE: If only a portion of the plate is used, for the 96-well plate format, we recommend reconstituting the oligos in 5 µL of water. Pipette 10-15 times and gently transfer 4.6 µL to the new plate where the RT will be performed. Do not add water to the Master Mix for the RT reaction (as indicated in step 2.1.5 below).

For the 384-well plate format, the volume for oligo reconstitution will depend on the sample volumes, and no additional water should be introduced into the master mix.

Before opening the seal, it is advisable to centrifuge the oligo plate, and **only the wells intended for use should be opened**. The seal must remain on the plate to prevent any potential cross-contamination.

2.1.1. Keep the RT plate on ice. Using a multichannel pipette, transfer the following volume of organoid lysate directly to the corresponding wells and pipette 3-5 times to ensure proper reconstitution of dried oligo-dT:

- **96WP:** 10 µL
- **384WP:** 5 µL

2.1.2. The appearance of red color in all wells indicates a proper and uniform reconstitution of oligos.

2.1.3. Carefully re-seal the RT plate and briefly spin it in the centrifuge.

2.1.4. Leave the RT plate on ice for 5 min.

2.1.5. Prepare the Master Mix for the RT reaction (+10%) as follows:

Reagent	96WP, µL		384WP, µL	
	Per well	96 wells +10%	Per well	384 wells +10%
RTB	5.0	528.0	2.5	1056
RTE	0.1	10.6	0.1	43
Water	4.9	517.5	2.4	1014
TOTAL	10.0	1056.1	5.0	2113

2.1.6. Keep the RT plate on ice and, using a multichannel pipette, transfer the following volume of Master Mix to each well containing the organoid lysate sample:

- **96WP**: 10 μ L
- **384WP**: 5 μ L

2.1.7. Carefully re-seal the RT plate and briefly spin it in the centrifuge.

2.1.8. Transfer the plate to the thermocycler and start **Program 1_RT**.

Safe stop: After this step, the RT plate can be kept at 4°C overnight.

2.2. Sample pooling and column purification

After pooling, the barcoded RT samples can be purified using either column-based Zymo Clean & Concentration Kit (Zymo, D4014) or SPRI magnetic beads (Beckman, A63881). Both approaches produce comparable outcomes and can be used interchangeably. Depending on the availability of 3rd party reagents and instruments, the corresponding method should be applied.

NOTE: The pool may contain some cell debris, which could block a column membrane during purification leading to a long waiting time. To avoid this, it is recommended to perform a pre-cleaning of the RT pool by passing it through the Zymo column (see below) before mixing it with 7x DNA Binding buffer.

The procedure of cDNA pre-cleaning and purification using the column-based method

After the cDNA from each well is pooled in a reservoir, mix it with a 7x volume of DNA binding buffer (Zymo, D4004-1-L). We strongly recommend using a vacuum manifold for the cDNA purification to avoid column membrane damage due to multiple centrifugation rounds. A high-capacity Zymo-Spin IIICG column (Zymo, C1006-50-G) is required to purify large volumes resulting from 384 sample pooling.

Plate format	Plasticware		Pooling volumes					
	Pipetting strategy	Zymo-Spin column type	First strand cDNA		DNA Binding Buffer		TOTAL	
			Per well, μ L	Per plate, mL	Per well, μ L	Per plate, mL	Per well, μ L	Per plate, mL
96WP	multichannel pipette or pipetting robot	I (#D4014)	20	1.92	140	13.44	160	15.36
384WP	pipetting robot	IIICG (C1006-50-G)	10	3.84	70	26.88	80	30.72

Table 1 Overview of the recommended pipetting strategy, plasticware, and reagent volumes to be used depending on the number of pooled samples

Preparation

- Make sure Zymo DNA Wash buffer has Ethanol added.

Procedure

2.2.1. According to **Table 1**, use a multichannel pipette, or pipetting robot to transfer the entire RT volume (20 μ L for **96WP**, 10 μ L for **384WP**) of each sample into a specific reservoir (25 mL or 100 mL).

2.2.2. Mix well the pool and transfer it to a falcon tube with a pipette.

2.2.3. For RT pool pre-cleaning, place a Zymo column in a new 2 mL tube, add 800 μ L of the collected pool, and briefly centrifuge.

2.2.4. Collect the cleaned flowthrough in a new falcon tube. Repeat step **2.2.3** until all the pool passes through the Zymo column. Discard the column.

2.2.5. Using a pipette, measure the volume of the pool after cleaning, transfer it to a 50 mL falcon tube and add 7x DNA Binding buffer accordingly (see **Table 1**). The color of the mix should turn yellow.

2.2.6. Connect the 25 mL funnel (Zymo, C1039-25) to a Zymo column suitable for purification volume (**Table 1**) and place it on a vacuum manifold.

- 2.2.7. Gently mix the cDNA in the binding buffer mixture and transfer it to 25 ml funnel using a pipetboy.
- 2.2.8. Turn on the vacuum pump and let the liquid pass through the column.
- 2.2.9. Transfer any remaining volume to the funnel. Do not let the membrane overdry.
- 2.2.10. After the entire pool mix has passed through the column, add 200 µL of DNA Wash buffer (with Ethanol added) directly to the membrane of the column.
- 2.2.11. Repeat step 2.2.10 once the wash buffer passed through the filter.
- 2.2.12. Remove the column from the vacuum manifold, put it in a 1.5 mL tube, and centrifuge for 1 min to remove leftovers from the washing buffer.
- 2.2.13. Depending on the Zymo-Spin column type used, perform the following:
 - For the type **I** column used with ≤96 samples (**96WP**), add 20 µL of water to the column matrix and incubate for 1 min.
 - For the type **IIICG** column used with 384 samples (**384WP**), add 38 µL of water to the column matrix and incubate for 1 min.
- 2.2.14. Transfer the column into a new labeled 1.5 mL tube and centrifuge for 30 sec.
- 2.2.15. Immediately proceed to step 2.3.

The procedure of cDNA purification using the SPRI bead-based method

Perform cDNA purification with SPRI magnetic beads using a 1:1 ratio of cDNA pool and beads slurry. The purification of large volumes (i.e., 2 mL of the pool from **96WP** and 4 mL from **384WP**) requires three to six 1.5 mL tubes and a corresponding magnetic stand (Permagen, MSR06).

If the volume of the pool is higher than 750 µL, split it equally in the required number of 1.5 mL tubes and add the identical volume of beads (i.e., a pool of 1 mL split in 2 tubes with 500 µl per tube and add 500µL of beads per tube).

NOTE: Please use SPRI-beads only if solution is clear and has no visible debris.

Preparation

Pre-warm beads at room temperature and vortex them vigorously before pipetting.

- 2.2.16. Pool the RT samples as described in Table 1.
- 2.2.17. Transfer the collected pool to a 2 mL or 15 mL tube, depending on the pooled volume. Consider that the final volume will be twice higher due to the addition of the beads.
- 2.2.18. Add pre-warmed beads in a 1:1 ratio (i.e., for 960 µL of pooled samples, add 960 µL of beads slurry), and mix by pipetting up and down ten times.
- 2.2.19. Incubate for 5 min at room temperature.
- 2.2.20. Place the tube on the magnetic stand, wait 5 min, and carefully remove and discard the supernatant.
- 2.2.21. To wash the beads, pipette 1 mL of freshly prepared 80% ethanol into the tube.
- 2.2.22. Incubate for 30 sec.
- 2.2.23. Carefully remove the ethanol without touching the bead pellet.
- 2.2.24. Repeat step 2.2.21 for a total of two washes.
- 2.2.25. Remove the tube from the magnetic stand and let the beads dry for 1-2 min.
- 2.2.26. Resuspend the beads in 37 µL of water and incubate for 1 min.
- 2.2.27. Place tubes on the magnetic stand, wait 5 min, and carefully transfer 35 µL of supernatant to a new tube to avoid bead carry-over.
- 2.2.28. Immediately proceed to step 2.3.

If the RT pool was split into several tubes at step 2.2.16, use one of the following options:

- **[Two tubes only]**, resuspend the beads in **both tubes** in 20 μL /tube, and combine both elutions in one tube;
- **[Two or more tubes]** resuspend the beads in the **first tube** in 40 μL of water. Keep other tubes closed to avoid over-drying of the beads. Transfer obtained elution to the next tube and resuspend beads. Repeat this step for every tube;
- **[Two or more tubes]**, resuspend **every** tube in 37 μL . Combine all elutions in one tube and perform one additional purification of the pool adding beads slurry accordingly to the pool volume (steps 2.2.18 - 2.2.27). Elute in 37 μL of water and collect 35 μL in a new tube.

2.3. Free primer digestion

It is recommended to perform non-incorporated primer digestion immediately after pooling.

Preparation

- Label 0.2 mL PCR tubes corresponding to the number of pools prepared.
- Thaw the **EXB** reagent at room temperature.
- Keep the **EXO** reagent on ice.
- Prepare program **2_FPD** on the thermocycler (set the lid at 90°C):

Step	Temperature, °C	Time
Incubation	37	30 min
Incubation	80	20 min
Keep	4	pause

Procedure

- 2.3.1. Depending on the cDNA volume, obtained from steps 2.2.15 or 2.2.27, transfer 17 μL or 35 μL of the eluate from each tube into a new labeled 0.2 mL PCR tube.
- 2.3.2. Prepare the EXO reaction mix as follows (with 10% excess):

Reagent	Per reaction (μL)	
	For 17 μL elution	For 35 μL elution
EXB	2.0	4.0
EXO	1.0	1.0
TOTAL	3.0	5.0

- 2.3.3. According to the table, transfer **3 μL** or **5 μL** of EXO reaction mix into each PCR tube with purified cDNA.
- 2.3.4. Mix by pipetting up and down 5 times.
- 2.3.5. Briefly spin down in the bench-top centrifuge.
- 2.3.6. Incubate in thermocycler Program **2_FPD**.
- 2.3.7. Proceed immediately to step 2.4 or keep the tube at 4°C overnight.

Safe stop: After this step, the tube(s) can be kept at 4°C overnight.

2.4. Second strand synthesis

At this step, double-stranded full-length cDNA is generated and purified using magnetic beads.

Preparation

- Pre-warm the SPRI beads at room temperature for ~30 min.
- Prepare 5 mL of 80% ethanol.
- Thaw the **SSB** reagent at room temperature and mix well before use.
- Keep the **SSE** reagent constantly on ice.
- Prepare program 3_SSS on the thermocycler (set the lid at 70°C):

Step	Temperature, °C	Time
Incubation	37	20 min
Incubation	65	30 min
Keep	4	pause

Procedure

- 2.4.1. Prepare the SSS reaction mix for the second strand synthesis as follows (with 10% excess):

Reagent	Per reaction (µL)	
	For 17 µL elution	For 35 µL elution
SSB	5.0	7.0
SSE	2.0	3.0
TOTAL	7.0	10.0

- 2.4.2. Transfer **7 µL** or **10 µL** of SSS reaction mix to the tube from step 2.3.7 and mix well by pipetting up and down 5 times.
- 2.4.3. Incubate in thermocycler Program 3_SSS.
- 2.4.4. Proceed immediately to step 2.4.5.

cDNA clean-up with SPRI beads

Perform the double-stranded cDNA purification with SPRI magnetic beads using a 0.6x ratio (i.e., 30 µL of bead slurry plus 50 µL of cDNA):

NOTE: Use pre-warmed beads and vortex them vigorously before pipetting.

- 2.4.5. Complement the final volume to 50 µL with water.
- 2.4.6. Add 30 µL of beads and mix by pipetting up and down 10 times.
- 2.4.7. Incubate for 5 min at room temperature.
- 2.4.8. Place the tube on the magnetic stand, wait 5 min, and carefully remove and discard the supernatant.
- 2.4.9. To wash the beads, pipette 200 µL of freshly prepared 80% ethanol into the tube.
- 2.4.10. Incubate for 30 sec.
- 2.4.11. Carefully remove the ethanol without touching the bead pellet.
- 2.4.12. Repeat step 2.4.9 for a total of two washes.
- 2.4.13. Remove the tube from the magnetic stand and let the beads dry for 1-2 min.
- 2.4.14. Resuspend the beads in 21 µL of water.
- 2.4.15. Incubate for 1 min.
- 2.4.16. Place the tube on the magnetic stand, wait 5 min, and carefully transfer 20 µL of supernatant into a new tube to avoid bead carry-over.
- 2.4.17. Use 2 µL to measure the concentration with Qubit.

Safe stop: At this step, the cDNA can be safely kept at -20°C for a few weeks.

2.5. Tagmentation

At this step, the full-length cDNA is tagmented using a Tn5 transposase pre-loaded with adapters for library amplification. It is recommended to use 20 ng of cDNA for tagmentation to obtain a higher library complexity with less PCR amplification.

Preparation

- Pre-warm the SPRI beads at room temperature for ~30 min.
- If needed, prepare fresh 5 mL of 80% ethanol.
- Thaw the **TAB** reagent at room temperature and mix well before use.
- Keep the **TE** reagent constantly on ice.
- Set the PCR machine at 55°C incubation.

Procedure

- 2.5.1. Prepare the Tagmentation mix on ice in a PCR tube. If several cDNA samples are tagmented, prepare a Master mix as follows (with 10% excess):

Reagent	Per reaction (µL), for cDNA inputs		
	≤ 9 ng	10-14 ng	15-20 ng
TAB	4.0	4.0	4.0
TE	1.0	2.0	3.0
cDNA	X µL	X µL	X µL
Water	15.0 – X	14.0 – X	13.0 – X
TOTAL	20.0	20.0	20.0

- 2.5.2. Keep the mix on ice and pipette up and down 10 times with the pipette set at 5 µL. Pay attention to mixing the reaction volume thoroughly.
- 2.5.3. Incubate for 7 min at 55°C in the PCR machine.
- 2.5.4. Immediately put the tube on ice and add 30 µL of water for the final volume of 50 µL.

Tagmented cDNA clean-up with SPRI beads

Purify the tagmented cDNA with SPRI magnetic beads using a 0.6x ratio.

NOTE: Use pre-warmed beads and vortex them vigorously before pipetting.

- 2.5.5. Add 30 µL of beads to 50 µL of tagmented cDNA and mix by pipetting up and down 10 times.
- 2.5.6. Incubate for 5 min at room temperature.
- 2.5.7. Place the tubes on the magnetic stand, wait 5 min, and carefully remove and discard the supernatant.
- 2.5.8. To wash the beads, pipette 200 µL of freshly prepared 80% ethanol into the tube.
- 2.5.9. Incubate for 30 sec.
- 2.5.10. Carefully remove the ethanol without touching the bead pellet.
- 2.5.11. Repeat steps 2.5.8 - 2.5.10 for a total of two washes.
- 2.5.12. Remove the tube from the magnetic stand and let the beads dry for 1-2 minutes.
- 2.5.13. Resuspend the beads in 21 µL of water.
- 2.5.14. Incubate for 1 min.
- 2.5.15. Place the tube on the magnetic stand, wait 5 min, and carefully transfer 20 µL of the supernatant into a new tube to avoid bead carry-over.
- 2.5.16. Proceed immediately to step 2.6.

2.6. Library indexing and amplification

At this step, 5' terminal fragments are amplified using the Unique Dual Indexing (UDI) adapter primers. The kit contains 4 Illumina-compatible primer pairs to generate the UDI libraries. The index sequences are indicated in [Table 4](#).

The number of amplification cycles required for library preparation typically ranges from 12 to 17. The precise number may depend on the RNA samples and the input cDNA amount used for tagmentation. To determine the optimal number of cycles, follow the Library quantification protocol below.

It is strongly recommended to perform the final library clean-up step twice to effectively remove primer dimer fragments.

Preparation

- Pre-warm the SPRI beads at room temperature for ~30 mins.
- Prepare 10 mL of 80% ethanol.
- Thaw the **LAM** reagent on ice and mix well before use.
- Thaw the required number of **MQ.UDI Adapters** at room temperature and briefly spin before use.
- Prepare the **Program 4_TN5AMP** (set the lid at 100°C) on the thermocycler (*The exact number of PCR cycles should be determined following the Library quantification protocol below)

Step	Temperature, °C	Time	Cycles
Initial denaturation	98	1 min	
Denaturation	98	10 sec	
Annealing	63	30 sec	15 or TBD*
Extension	72	1 min	
Final extension	72	3 min	
Keep	4	pause	

Procedure

2.6.1. Prepare the PCR amplification reaction as follow:

Reagent	Per reaction (µL)
LAM	25.0
MQ.UDI Adapter	5.0
Tagmented cDNA	20.0
TOTAL	50.0

2.6.2. Pipette up and down 5 times.

2.6.3. Put the tube in the PCR machine and do one of the following:

- Determine the optimal number of amplification cycles using real-time PCR (highly recommended; follow steps [2.6.5](#) - [2.6.12](#)); or
- start **Program 4_TN5AMP** with the default 15 PCR cycles (less preferable).

2.6.4. After the PCR is done, proceed to the library bead clean-up ([2.6.13](#)).

2.6.5. In the case of the cycle number optimization, perform 5 cycles of library preamplification (**Program 4_TN5AMP**).

2.6.6. Place the tube on ice.

2.6.7. Use a 5 µL aliquot from the PCR reaction to prepare a qPCR reaction mix in the appropriate PCR tube or plate as follows:


Reagent	Per reaction, µL
5 cycles pre-amplified library	5.0
SYBR 100x*	0.1
LAM	2.5
Water	2.4
TOTAL	10.0

2.6.8. Place the tube in the qPCR machine and run with the following program:

Step	Temperature	Time	Cycles
Initial denaturation	98°C	1 min	
Denaturation	98°C	10 sec	
Annealing	63°C	30 sec	25
Extension + acquisition	72°C	1 min	

2.6.9. Determine the cycle number depending on the growth curve in the multicomponent plot, as in Figure 2.

Figure 2 Determination of the additional number of amplification cycles with qPCR



2.6.10. Place the tube containing the remaining 45 μ L of the pre-amplified library in the PCR machine.

2.6.11. Set the number of amplification cycles determined in step 2.6.9.

2.6.12. Typically, 10 ng of cDNA is sufficient to obtain 20-40 ng of DNA library after 13-14 cycles of amplification. Table 2 shows the approximate amount of amplification cycles and the expected library yield.

Tagmented cDNA, ng	# PCR cycles	Expected library yield, ng
5	14-15	
10	13-14	20-40
20	11-12	

Table 2 Expected yield of BRB-seq libraries from different amounts of tagmented cDNA

Indexed cDNA library clean-up with SPRI beads

Purify the final cDNA library with SPRI magnetic beads using a 0.7x ratio (35 μ L of bead slurry for 50 μ L cDNA library).

NOTE: Use pre-warmed beads and vortex them vigorously before pipetting.

2.6.13. Adjust the library volume to 50 μ L with water.

2.6.14. Add 35 μ L of beads and mix by pipetting up and down 10 times.

2.6.15. Incubate for 5 min at room temperature.

2.6.16. Place the tubes on the magnetic stand, wait 5 min, carefully remove, and discard the supernatant.

2.6.17. To wash the beads, pipette 200 μ L of freshly prepared 80% ethanol into the tube.

2.6.18. Incubate for 30 sec.

2.6.19. Carefully remove the ethanol without touching the bead pellet.

2.6.20. Repeat step 2.6.17 for a total of two washes.

2.6.21. Remove tubes from the magnetic stand and let the beads dry for 1-2 min.

2.6.22. Resuspend the beads in 21 μ L of water.

2.6.23. Incubate for 1 min.

2.6.24. Place tubes on the magnetic stand, wait 5 minutes, and carefully remove 20 μ L of supernatant into a new tube to avoid bead carry-over.

2.6.25. Perform the bead clean-up once again by repeating the procedure from step 2.6.13.

Safe stop: At this step, the cDNA libraries can be safely kept at -20°C for a few weeks.

2.7. Library quality control

Pooled library quality control

Before sequencing, the libraries should be subjected to fragment analysis (with Fragment analyzer, Bioanalyzer, or TapeStation) and quantification (with Qubit). This information is required to assess the molarity of the libraries and prepare the appropriate library dilution for sequencing. A successful library contains fragments in the range of 300 – 1000 bp with a peak at 400-700 bp; see Figure 3 for an example of a standard BRB-seq library profile.

Importantly, libraries with primer dimer peaks at 180 bp and ranging at 250 – 290 bp will likely produce lower-quality sequencing data with a reduced proportion of mappable reads (Figure 4). Therefore, it is strongly recommended to remove those peaks by performing an additional round of SPRI beads purification with the 0.7x ratio (see step 2.6.13).

Undertagged libraries have a broader fragment range distribution with a peak at 1000-1200 bp (Figure 5). Only a fraction of such libraries contains fragments that can be efficiently sequenced; therefore, it is recommended to re-tagment the cDNA for the best results.

Library quantification can also be done unbiasedly by qPCR using standard Illumina library quantification kits (i.e., KAPA HiFi, Roche).

Pre-sequencing library QC:

- Use 2 μ L of the library to measure the concentration with Qubit.
- Use 2 μ L of the library to assess the profile with the Fragment Analyzer instrument or similar.
- If necessary, re-purify the libraries by following steps 2.6.13 – 2.6.24 to remove the peaks <300 bp.

Figure 3 A successful library profile with fragments between 300-1000 bp




Figure 4 An example of an over-tagmented library profile with a peak at 290 bp and an adapter peak at 160 bp





Figure 5 An example of an under-tagmented library profile with a major peak > 1000 bp



Assessing uniformity of read distribution across the samples

For projects involving highly heterogeneous RNA samples, it is recommended to validate the uniformity of read coverage across the samples by shallow library sequencing (see step 2.2). This approach ensures that every sample will obtain enough reads required for the analysis. BRB-seq libraries can be added as spike-ins to the compatible sequencing run (see Part 3). For this validation 0.5-1M sequencing reads per library is sufficient to assess the fraction of reads attributed to every sample.

Part 3. LIBRARY SEQUENCING

The libraries prepared with the MERCURIUS™ BRB-seq kit carry Illumina- and AVITI-compatible adapter sequences. They can be processed on any Illumina instrument (e.g., HiSeq, NextSeq, MiSeq, iSeq, and NovaSeq) or in the Element AVITI System with Adept Workflow.

The MERCURIUS™ BRB-seq libraries are Unique Dual-Indexed and can potentially be pooled in a sequencing run with other libraries if the sequencing structure is compatible. Please refer to [Table 3](#) for the optimal sequencing structure and [Table 4](#) for the list of i5 and i7 index sequences.

Given the BRB-seq library structure, the optimal number of cycles for Read 1 is 28 (and 29 for AVITI). The following cycles, 29-60, will cover the homopolymer sequence, which may result in a significant drop of Q30 values reflecting sequencing quality. Therefore, using standard Illumina or AVITI runs setups (e.g., 100 PE or 150 PE) is not recommended.

Read	Length (cycles)		Comment
	for Illumina	for AVITI	
Read 1	28	29	Sample barcode (14 nt) and UMI (14 nt); +1 extra base for AVITI
Index 1 (i7) read	8	8	Library Index
Index 2 (i5) read	8*	8*	Library Index (*optional and valid for UDI libraries)
Read 2	60-90	101	Gene fragment

Table 3 Sequencing structure of BRB-seq libraries

The Unique Dual Indexing (UDI) strategy ensures the highest library sequencing and demultiplexing accuracy and complies with the best practices for Illumina sequencing platforms. UD-indexed libraries have distinct index adapters for i7 and i5 index reads ([Table 4](#)).

Name	Type	i7 index sequence	i5 index sequence Forward Workflow	i5 index sequence Reverse Workflow
MQ.UDI.1	UDI (i7/i5)	TAAGGCCGA	TATAGCCT	AGGCTATA
MQ.UDI.2	UDI (i7/i5)	CGTACTAG	ATAGAGGC	GCCTCTAT
MQ.UDI.3	UDI (i7/i5)	AGGCAGAA	CCTATCCT	AGGATAGG
MQ.UDI.4	UDI (i7/i5)	GCGTTGGA	TTGGACTT	AAGTCCAA

Table 4 UDI adapter sequences

NOTE: Sequencing depth

1. The recommended sequencing depth is 1-5 Mio reads per sample. Deeper sequencing can also be performed to detect very lowly expressed genes.
2. If only one library is sequenced in a flow cell, the Index reads can be skipped.
3. The loading molarity for the library depends on the type of sequencing instrument (see 3.1 and 3.2) and should be discussed with the sequencing facility or an experienced person.

3.1. Sequencing with Illumina instruments

The loading concentration for the Illumina instruments is indicated in [Table 5](#). Please refer to [Appendix 2](#) for the list of Illumina instruments with forward or reverse workflow.

Instrument	Final loading concentration	PhiX
MiSeq	20 pM	1 %
iSeq	100 pM	1 %
NextSeq 500/550/550Dx	2.2 pM	1 %
NextSeq 2000, manual denature	85 pM	1%
NextSeq 2000, onboard denature	850 pM	1%
NovaSeq Standard Workflow*	160 pM	1 %
NovaSeq XP Workflow	100 pM	1 %
HiSeq4000	270 pM	1 %

* - adjusted molarity for BRB-libraries sequencing. We recommend a prior dilution of the libraries to 0.8 nM before denaturation.

Table 5 Reference loading concentrations for various Illumina instruments

3.2. Sequencing with Element AVITI instrument

For the most optimal results, the MERCURIUS™ BRB-seq libraries can be sequenced with the Element Biosciences AVITI System using Cloudbreak AVITI 2x75 High Output sequencing kits (#860-00004).

Libraries must be converted with the Adept PCR-Plus module (#830-00018) for linear loading (Table 6).

Type	Loading molarity, pM	Library starting amount for denaturation, nM	PhiX control	PhiX, %
Cloudbreak	14	1*	PhiX Control Library, Adept	2 %

* - requires 4nM of library before conversion

Table 6 Loading concentration for Cloudbreak AVITI 2x75 High Output sequencing kit.

NOTE: Sequencing depth

Please note that the **Cloudbreak AVITI 2x75 High Output sequencing** yields 1 Mio reads. Therefore, for the 384-sample library, the sequence depths would be limited to a maximum of 2.5 Mio reads/sample.

Part 4. SEQUENCING DATA PROCESSING

Following Illumina sequencing and standard library index demultiplexing, the user obtains raw read1 and read2 *fastq* sequencing files (e.g., *mylibrary_R1.fastq.gz* and *mylibrary_R2.fastq.gz*).

This section explains how to generate ready-for-analysis gene and UMI read count matrices from raw *fastq* files.

This section explains how to generate ready-for-analysis gene, and UMI read count matrices from raw *fastq* files.

To obtain the data ready for analysis, the user needs to align the sequencing reads to the genome and perform the gene/UMI read count generation, which can be done in parallel with the sample demultiplexing.

For manual data processing, the user requires a terminal and a server or powerful computer with an installed set of standard bioinformatic tools.

4.1. Required software

- [fastQC](#) (version v0.11.9 or greater). Software for QC of *fastq* or *bam* files. This software is used to assess the quality of the sequencing reads, such as the number of duplicates, adapter contamination, repetitive sequence contamination, and GC content. The software is freely available from <https://www.bioinformatics.babraham.ac.uk/projects/fastqc/>. The website also contains informative examples of [good](#) and poor-quality data.
- [STAR](#)solo from STAR (version 2.7.9a). Software for read alignment on reference genome (Dobin et al., 2013). It can be downloaded from [Github](https://github.com/alexdobin/STAR) (<https://github.com/alexdobin/STAR>). STAR can only be run on UNIX systems and requires:
 - x86-64 compatible processors
 - 64-bit Linux or Mac OS X.
 - ~30-40Gb of RAM
- [FastReadCounter](#) (v.1.1 or greater). Software for counting genome-aligned reads for genomic features. github.com/DeplanckeLab/FastReadCounter
- [Picard](#) (v.2.17.8 or greater) and [Samtools](#) (v.1.9 or greater). Collections of command-line utilities to manipulate with BAM files. Note: Picard requires [Java version 8 or higher](#) to be installed.
- R Software (version 3 or greater).
- (Optional) [BRBseqTools](#) (version 1.6). The software suite for processing BRB-seq libraries is available at <https://github.com/DeplanckeLab/BRB-seqTools>.

4.2. Data processing

4.2.1. Merging *fastq* files from individual lanes and/or libraries (Optional)

- 4.2.1.1 Depending on the type of instrument used for sequencing, one or multiple R1/R2 *fastq* files per library may result from individual lanes of a flow cell. The *fastq* files from individual lanes should be merged into single *R1.fastq* and single *R2.fastq* files to simplify the following steps. This is an example of *fastq* files obtained from HiSeq 4 lane sequencing:

```
> mylibrary_L001_R1.fastq.gz, mylibrary_L002_R1.fastq.gz,
  mylibrary_L003_R1.fastq.gz, mylibrary_L004_R1.fastq.gz
> mylibrary_L001_R2.fastq.gz, mylibrary_L002_R2.fastq.gz,
  mylibrary_L003_R2.fastq.gz, mylibrary_L004_R2.fastq.gz
```

- 4.2.1.2 To merge the *fastq* files from different lanes use a `cat` command in a terminal. This will generate two files: *mylibrary_R1.fastq.gz* and *mylibrary_R2.fastq.gz*, containing the information of the entire library.

```
> cat mylibrary_L001_R1.fastq.gz mylibrary_L002_R1.fastq.gz
  mylibrary_L003_R1.fastq.gz mylibrary_L004_R1.fastq.gz >
mylibrary_R1.fastq.gz
> cat mylibrary_L001_R2.fastq.gz mylibrary_L002_R2.fastq.gz
  mylibrary_L003_R2.fastq.gz mylibrary_L004_R2.fastq.gz >
mylibrary_R2.fastq.gz
```

- 4.2.1.3 Move these 2 *fastq* files into a new folder, which will be referenced in this manual as `$fastqfolder`.

NOTE: This step can also be done if you sequenced your library in multiple sequencing runs.

Warning: The order of merging files should be kept the same (e.g., L001, L002, L003, L004, not L002, L001 ...) to avoid issues when demultiplexing the samples.

4.2.2. Sequencing data quality check

4.2.2.1 Run fastQC on both R1 and R2 fastq files. Use `--outdir` option to indicate the path to the output directory. This directory will contain HTML reports produced by the software.

```
> fastqc --outdir $QCdir/ mylibrary_R1.fastq.gz
> fastqc --outdir $QCdir/ mylibrary_R2.fastq.gz
```

4.2.2.2 Check fastQC reports to assess the quality of the samples (see Software and materials).

NOTES:

- The report for the R1 *fastq* file may contain some "red flags" because it contains barcodes/UMIs. Still, it can provide useful information on the sequencing quality of the barcodes/UMIs.
- The main point of this step is to check the R2 *fastq* report. Of note, *per base sequence content* and *kmer content* are rarely green. If some *adapter contamination* or *overrepresented sequence* is detected in the data, it may not be an issue (if the effect is limited to <10~20%). These are lost reads, but most of them will be filtered out during the next step.

4.2.3. Preparing the reference genome

The *fastq* files must be aligned (or "mapped") on a reference genome. The [STAR](#) (Dobin et al., 2013¹) aligner is one of the most efficient tools for RNA-seq reads mapping. It contains a "soft-clipping" tool that automatically cuts the beginning or the end of reads to improve the mapping efficiency, thus allowing the user to skip the step of trimming the reads for adapter contamination. Moreover, STAR has a mode called STARsolo, designed to align multiplexed data (such as BRB-seq) and directly generate count matrices.

The STAR aligner requires a genome assembly together with a genome index file. The index file generation is a time-consuming process that is only performed once on a given genome assembly so that it can be completed in advance and the index files can be stored on the server for subsequent analyses.

4.2.3.1 Download the correct genome assembly fasta file (e.g. Homo_sapiens.GRCh38.dna.primary_assembly.fa) and gene annotation file in gtf format (e.g. Homo_sapiens.GRCh38.108.gtf) from Ensembl or UCSC repository. Below is an example for a human assembly:

```
> wget https://ftp.ensembl.org/pub/release-108/fasta/homo_sapiens/dna/Homo_sapiens.GRCh38.dna.primary_assembly.fa.gz
> gzip -d Homo_sapiens.GRCh38.dna.primary_assembly.fa.gz # unzip
> wget https://ftp.ensembl.org/pub/release-108/gtf/homo_sapiens/Homo_sapiens.GRCh38.108.gtf.gz
> gzip -d Homo_sapiens.GRCh38.108.gtf.gz # unzip
```

NOTE: It's recommended to download the primary_assembly fasta file when possible (without the 'sm' or 'rm' tags). If not available, download the top_level assembly. For the gtf, download the one that does not have the 'chr' or 'abinitio' tags.

4.2.3.2 Use STAR to create an index for the genome assembly. Indicate the output folder name containing the index files using `--genomeDir` option:

```
> STAR --runMode genomeGenerate --genomeDir /path/to/genomeDir --
genomeFastaFiles Homo_sapiens.GRCh38.dna.primary_assembly.fa --sjdbGTFfile
Homo_sapiens.GRCh38.108.gtf --runThreadN 8
```

NOTES:

¹ Alexander Dobin, Carrie A. Davis, Felix Schlesinger, Jorg Drenkow, Chris Zaleski, Sonali Jha, Philippe Batut, Mark Chaisson, Thomas R. Gingeras, STAR: ultrafast universal RNA-seq aligner, *Bioinformatics*, Volume 29, Issue 1, January 2013, Pages 15–21, <https://doi.org/10.1093/bioinformatics/bts635>

- The `--runThreadN` parameter can be modified depending on the number of cores available on your machine. The larger this number is, the more parallelized/fast the indexing will be.
- Depending on the genome assembly, STAR can use up to 32-40Gb of RAM. So, you should use a machine that has this RAM capacity.

4.2.4. Aligning to the reference genome and generation of count matrices

After the genome index is created, both R1 and R2 *fastq* files can be aligned to this reference genome. For this step, use the “solo” mode of STAR, which not only aligns the reads to the reference genome but also creates the gene read count and UMI (unique molecular identifier) count matrices.

The following parameters should be adjusted according to the sequencing information:

- `--soloCBwhitelist`: a text file with the list of barcodes (one barcode sequence per lane) which is used by STAR for demultiplexing. This file is provided according to version of the MERCURIUS kit used. Example of “[barcodes_96_V5D_star.txt](#)”:

```
> TACGTTATTCCGAA
> AACAGGATAACTCC
> ACTCAGGCACCTCC
> ACGAGCAGATGCAG
```

- `--soloCBstart`: Start position of the barcode in the R1 *fastq* file, equal to 1.
- `--soloCBlen`: Length of the barcode. This value should match the length of the barcode sequence in the file specified by `--soloCBwhitelist`. The length of the barcode depends on the version of the oligo-dT barcodes provided in the kit. For the barcode plate set V5, the default value is 14.
- `--soloUMIstart`: Start position of the UMI, it's `soloCBlen + 1` since the UMI starts right after the barcode sequence.
- `--soloUMIlen`: The length of UMI. This parameter depends on the version of the oligo-dT barcodes in the kit and the number of sequencing cycles performed for Read1. For the barcode plate set V5 the default value is 14.
- `--readFilesIn`: name and path to the input *fastq* files.

The order of the *fastq* files provided in the script is important. The first *fastq* must contain genomic information, while the second the barcode and UMI content. Thus, files should be provided for STARsolo in the following order: `--readFilesIn mylibrary_R2 mylibrary_R1`.

- `--genomeDir`: a path to the genome indices directory generated before (`$genomeDir`).

Output count matrix parameters:

By default, STARsolo produces a UMI count matrix, i.e., containing unique non-duplicated reads per sample for each gene. This type of count data is a standard for single-cell RNA-seq analysis. For bulk RNA-seq analysis, a gene read count matrix is usually used. The following parameters will enable the generation of the output of interest.

`--soloUMI dedup NoDedup`, will generate a read count matrix output

`--soloUMI dedup NoDedup 1MM_Directional`, will generate both UMI and read count matrices in *mtx* format.

This step will output *bam* files and count matrices in the folder `$bamdir`.

```
> STAR --runMode alignReads --outSAMmapqUnique 60 --runThreadN 8 --
outSAMunmapped Within --soloStrand Forward --quantMode GeneCounts --
outBAMsortingThreadN 8 --genomeDir $genomeDir --soloType CB_UMI_Simple --
soloCBstart 1 --soloCBlen 14 --soloUMIstart 15 --soloUMIlen 14 --
soloUMI dedup NoDedup 1MM_Directional --soloCellFilter None --soloCBwhitelist
barcodes.txt --soloBarcodeReadLength 0 --soloFeatures Gene --
outSAMattributes NH HI nM AS CR UR CB UB GX GN sS sQ sM --
outFilterMultimapNmax 1 --readFilesCommand zcat --outSAMtype BAM
SortedByCoordinate --outFileNamePrefix $bamdir --readFilesIn
mylibrary_R2.fastq.gz mylibrary_R1.fastq.gz
```

The demultiplexing statistics can be found in the “*bamdir/Solo.out/Barcodes.stats*” file.

The alignment quality and performance metrics can be found in the “*bamdir/Log.final.out*” file.

NOTE: The most important statistic at this step is the proportion of “Uniquely mapped reads” which is expected to be greater than 70% (for human, mouse or drosophila).

4.2.5. Generating the count matrix from .mtx file

STARsolo will generate a count matrix (*matrix.mtx* file) located in the *bamdir/Solo.out/Gene/raw* folder. This file is a sparse matrix format that can be transformed into a standard count matrix using an R script provided below:

```
> #Myscript.R
> library(data.table)
> library(Matrix)
> matrix_dir <- "$bamdir/Solo.out/Gene/raw"
> f <- file(paste0(matrix_dir, "matrix.mtx"), "r")
> mat <- as.data.frame(as.matrix(readMM(f)))
> close(f)
> feature.names = fread(paste0(matrix_dir, "features.tsv"), header = FALSE,
stringsAsFactors = FALSE, data.table = F)
> barcode.names = fread(paste0(matrix_dir, "barcodes.tsv"), header = FALSE,
stringsAsFactors = FALSE, data.table = F)
> colnames(mat) <- barcode.names$V1
> rownames(mat) <- feature.names$V1
> fwrite(mat, file = umi.counts.txt, sep = "\t", quote = F, row.names = T,
col.names = T)
```

The resulting UMI/gene count matrix can be used for a standard expression analysis following conventional bioinformatic tools.

4.2.6. Generating the read count matrix with per-sample stats (Optional)

Given a multiplex BAM file obtained with STARsolo and a set of barcodes, the software FastReadCounter produces a read count matrix with per-sample statistics with the following code:

```
> #!/bin/bash
>
> gtf_file=homo_sapience.gtf          ### GTF genome annotation file
> output_folder=counts/              ### Name of the final count output file
> bam_dir=myspath/bam_demult         ### Directory with demultiplexed bam
files
> barcode_file=V5D_96_frc.txt        ### Barcode reference file
>
> FastReadCounter-1.0.jar" --bam ${bam_dir}/${bam_dir}.bam \
>                               --gtf "${gtf_file}" \
>                               --umi-dedup none \
>                               --barcodeFile ${barcode_file} \
>                               -o ${output_folder}
```

The resulting read count matrices can be used for subsequent gene expression analysis using established pipelines and tools.

NOTE: Please contact us at info@alitheagenomics.com in case you don't have the barcode sequences (in your email, please indicate the name of the barcode set and the PN of the barcode module)

4.2.7. Demultiplexing bam files (Optional)

Generation of demultiplexed bam files, i.e., individual bam files for each sample, might be needed in some cases, for example, for submitting the raw data to an online repository that does not accept multiplexed data (for example, GEO or ArrayExpress), or for running an established bulk RNA-seq data analysis pipeline.

For this purpose, the Picard tool can be used with the following parameters:

- `$out_dir`, The output directory for demultiplexed bam files
- `$path_to_bam`, the path to multiplexed single bam file
- `$barcode_brb.txt`, the tab-delimited file containing 2 columns: `sample_id` and `barcode seq`. Example of [barcode_96_V5D_brb.txt](#):

```
> Sample1      TACGTTATTCCGAA
> Sample2      AACAGGATAACTCC
> Sample3      ACTCAGGCACCTCC
```



```
> Sample4      ACGAGCAGATGCAG
```

NOTE: This file differs from the list of barcode files provided to STAR.

Run the following Picard script

```
> #!/bin/bash
> demultiplexed_bam_out_dir=$out_dir
> input_bam=$path_to_bam
> barcode_info=$barcode_brb.txt
>
> while IFS=$'\t' read -r -a line
> do
>     sample_id="${line[0]}"
>     tag_value="${line[1]}"
>
>     java -jar /path/to/picard.jar FilterSamReads \
>         I=${input_bam} \
>         O=${demultiplexed_bam_out_dir}/${sample_id}.bam \
>         TAG=CR TAG_VALUE=${tag_value} \
>         FILTER=includeTagValues
> done < "$barcode_info"
```

NOTE: Please contact us at info@alitheagenomics.com in case you don't have the barcode sequences (in your email, please indicate the name of the barcode set and the PN of the barcode module)

Appendix 1. ERCC Spike-In Control

The current protocol includes the addition of External RNA Controls Consortium (ERCC) Spike-Ins to the lysate buffer.

Prepare a 1:100 dilution of the ERCC RNA Spike-In mix in nuclease-free water. Mix 990 μ L of pre-chilled water with 10 μ L of ERCC. Pipette well and aliquot the dilution into 50 μ L aliquots, keeping them at -20°C.

The working solution of 1x Organoid Lysis Buffer with ERCC Spike-In controls consists of the following:

Reagent	96WP, μ L		384WP, μ L	
	Per well	96 wells +10%	Per well	384 wells +10%
DIR	15	1575	5	2200
OLB	10	1050	3.4	1496
INH	2.5	263	0.8	352
Water	1.3	136	0.7	308
ERCC* (1:100)	1.2	126	0.1	44
TOTAL	30	3150	10	4440

*The final ERCC is 1:1000 in a 384-type well (equal to 50 ng of RNA/well) and 1:250 in a 96-type well (150-200 ng of RNA/well)

Organoid Lysis Buffer (OLB) preparation with ERCC

1. Thaw the DIR, OLB, and ERCC tubes on ice and avoid their long-term storage.
2. Keep the nuclease-free water on ice to maintain a cold temperature.
3. Spin down all the tubes before pipetting.
4. First, add the water to a 15-mL falcon tube, then the DIR, OLB, OIN, and the ERCC (in this particular order).
5. Pipette the prepared mix a few times and briefly spin the tube. Keep it on ice until further use.
6. Follow the main protocol for the organoid lysis procedure ([step 1.4.1](#))

Appendix 2. Compatible Illumina instruments

Illumina instruments can use two workflows for sequencing i5 index (see the details in [Indexed Sequencing Overview Guide](#) on Illumina's website).

Forward strand workflow instruments:

- NovaSeq 6000 with v1.0 reagents
- MiSeq with Rapid reagents
- HiSeq 2500, HiSeq 2000

Reverse strand workflow instruments:

- NovaSeq 6000 with v1.5 reagents
- iSeq 100
- MiniSeq with Standard reagents
- NextSeq
- HiSeq X, HiSeq 4000, HiSeq 3000

Contact information

Alithea Genomics SA

Route de la Corniche 8
1066 Épalinges
Switzerland
Tel: +41 78 830 31 39

www.alitheagenomics.com
info@alitheagenomics.com